

Trifluoroacetic Acid

CAS No. 76-05-1

INTRODUCTION

Since it was discovered in 1922, trifluoroacetic acid has proved to be a significant chemical with very distinctive properties. While it smells like acetic acid, it is a much stronger acid - much more reactive and a better solvent - with an interesting list of uses. Trifluoroacetic acid (TFA) provides a way to introduce trifluoromethyl groups into more complex molecules. This compound has been used in the production of pharmaceutical and agricultural chemicals, as well as in many other specialized applications. It is used in peptide synthesis and as a solvent and catalyst in polymerization and condensation reactions. Halocarbon is proud to have been the first to manufacture this compound and remains the leading producer.

PROPERTIES


CAS Number	76-05-1
Molecular Formula	CF ₃ COOH
Molecular Weight	114.02
Freezing Point, °C	-15.0
Boiling Point, °C	72
Density, g/mL, 25°C	1.478
Refractive Index, <i>n</i> (25/d)	1.283
Heat of Vaporization, kJ/mol (Btu/lb)	33.5 (126)

¹ For more detailed information, see "The Chemistry of Nonaqueous Solvents" edited by J.J. Lagowski, Volume VB, Academic Press, New York, 1978, chapter on Trifluoroacetic Acid by John B. Milne, pp 1-52.

² Water Azeotrope with TFA (20.6 wt. % H₂O) is 105.5 °C

CHEMISTRY

Trifluoroacetic acid (TFA) reacts with active hydrogen (HZ, such as hydroxy and amino groups):




When using TFA, water is the coproduct, which works well in aqueous systems or systems tolerant of water.

The reactions of TFA are centered around the carbonyl group. The trifluoromethyl group, on the other hand, is one of the most stable organic structures. It is inert toward practically all oxidizing, reducing and hydrolyzing conditions and in vivo metabolism. The carbonyl group provides a way to attach the trifluoromethyl group to other molecules.

The usual carboxylic acid derivatives such as esters and amides can be made readily. They differ from the corresponding acetic acid derivatives in that the esters hydrolyze more readily and the amides of primary and secondary amines are sufficiently acidic to dissolve in 5% aqueous alkali.

TFA adds to olefins and acetylenes. The saturated ester from the olefin addition can be hydrolyzed readily to the corresponding alcohol. In fact, TFA can be used to promote hydration of olefins which are too sensitive to be hydrated in the usual sulfuric acid procedure. The vinyl ester made from acetylene polymerizes and the polymer is easily hydrolyzed to polyvinyl alcohol. Some of these general reactions are shown in the scheme below.




TFA is used as an esterification and transesterification promoter. In the presence of another carboxylic acid, that acid will be preferentially esterified. The mixed anhydride can be made from the other acid anhydride and TFA.

TFA can be used to promote the acylation of aromatic and unsaturated compounds to form the corresponding ketones and a variety of acyl-containing compounds.

TFA can be used to protect active amino groups temporarily. Esters of amino acids protected this way have sufficient volatility to be distilled or analyzed by GLC. Removal of the TFA is easily accomplished without changing chirality or other active groups in the molecule.

In some automated peptide reactions, TFA is added to deprotect amino groups sequentially as the molecule is being built. For this use, the solvent characteristics, acidity and volatility of TFA play a role. Our highly purified, low-

residue BioGrade™ TFA helps assure purity in the final protein. An example where the “BOC”-protecting group is removed from an amine is shown below. (WO03048158, 2003, Bristol-Myers Squibb)



Making trifluoroacetyl derivatives permits the use of ^{19}F NMR. Analysis of the chemical shifts can aid in determining the structure of the original molecule.

TFA can be mixed with other acids, or Lewis acids, to generate modified catalysts useful in olefin/alkane reactions in petroleum refining steps. These catalysts reduce byproduct sludges and increase yields when compared to the more conventional catalysts.

Peroxytrifluoroacetic acid is readily made from TFA and hydrogen peroxide. The peroxy acid is an oxidizer of amines, olefins, ketones, oximes and aromatics. The combination of the peroxy acid and boron trifluoride is an excellent source of positive hydroxyl groups for the hydroxylation of aromatic compounds.

The peroxy acid will dissolve metals not normally dissolved by mineral acids (Ag, Bi, Cu, Hg, In, Pb, Tl, etc.). This is useful for analytical purposes. Another analytical use for TFA is as the solvent during oxidations of metal ions by chromic acid, permanganate and chlorine. TFA is both inert and a good solvent in these cases.

TFA is miscible with most organic solvents and it will dissolve limited quantities of lower alkanes. Fluorinated solvents, even perfluorocarbons, are completely miscible with TFA. Polyamides and polyesters are soluble and the solutions can be used for analytical purposes or to perform reactions on the polymers. The solutions may also be used to apply coatings of the polymer or to spin fibers at low temperatures.

HANDLING AND SAFETY

Materials of Construction

TFA with less than 0.1% water is not corrosive to aluminum, nickel, monel, steels and higher alloys (see table). It can also be handled in glass and plastics such as polyethylene, polypropylene, PTFE, PCTFE and others, but it reacts with oxygen-containing plastics. Typically, TFA is shipped in polyethylene carboys, lined drums or larger stainless steel containers.

CORROSION by TFA at 71°C (30-day test)									
Metal	Nickel	Monel	SS304	SS316	Carp.20	Incoloy 825	Inconel	Aluminum	Copper
Mils/yr	0.4	0.7	1	0.07	0.06	0.1	0.4	0.03	5

Safety

TFA will cause severe burns to the skin, eyes, mucous membranes and other exposed tissues because of its strong acidity and fast tissue penetration. Special precautions and protective clothing, including rubber gloves and face masks, should be worn when handling it. The best treatment for spills on the skin is immediate and continuous flooding with water on the affected area.

TFA does not exhibit poisoning due to fluoride ion or any toxicity resembling that of monofluoroacetic acid.

For additional safety information refer to our MSDs. The information given above is intended for general reference only.